

INT.L DISPOSABLE NEONATAL ELECTRODES CATALOGUE

Valid from 01/03/2024

Code	Item Description	Use	Qty per pack	Qty per carton	
------	------------------	-----	--------------	----------------	--

ECG disposable neonatal electrodes

62.100.14	20mm disposable surface electrode with 35cm leadwire (B/R/Y) and 4mm socket (3 electrodes per pouch)	Neonatal	120	120		
62.100.01	24mm disposable surface electrode with 60cm leadwire (B/R/Y) and 1,5mm touch proof connector (3electrodes per pouch)	Neonatal	120	120		
62.100.02	24mm disposable surface electrode with 60cm leadwire (B/R/Y) and 4mm socket (3 electrodes per pouch)	Neonatal	120	120		
62.100.03	24mm disposable surface electrode with 90cm leadwire (B/R/Y) and 1,5mm touch proof connector (3 electrodes per pouch)	Neonatal	120	120		
62.100.05	20mm disposable surface electrode with 60cm leadwire (B/R/Y) and 1,5mm touch proof connector (3 electrodes per pouch)	Neonatal	120	120		
62.100.07	20mm disposable surface electrode with 90cm leadwire (B/R/Y) and 1,5mm touch proof connector (3 electrodes per pouch)	Neonatal	120	120		
62.111.01	24mm disposable surface electrode with 60cm leadwire (radiolucent) and 1,5mm touch proof connector (3 electrodes per pouch)	Neonatal	120	120		
62.111.03	20mm disposable surface electrode with 60cm leadwire (radiolucent) and 4mm socket (3 electrodes per pouch)	Neonatal	120	120		
64.140.05	24mm disposable surface electrode with snap connection (500 electrodes/carton - 25 electrodes/pouch)	Neonatal	500	500		
64.140.07	35x45mm disposable surface electrode with snap connection (80 electrodes/carton - 4 electrodes/pouch)	Neonatal	500	500		

The present catalogue nullifies and substitutes all the previous one.

APPROVED
21/02/2024